

Product Datasheet

Human Transmembrane Protein 119 ELISA Kit (orb1212826)

Description	Human Transmembrane Protein 119, TMEM119 ELISA Kit
Reactivity	Human
Range	20-7000ng/L
Tested Applications	ELISA
Target	TMEM119
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,203
Uniprot ID	Q4V9L6
Solubility (25°C)	8000ng/L
Sensitivity	9.25ng/L
Entrez	338773
Expiration Date	Please enquire.