

Product Datasheet

Human Transcriptional Repressor Ctcf ELISA Kit (orb1212830)

Description	Human Transcriptional Repressor Ctcf, CTCF ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	CTCF
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	82,785
Uniprot ID	P49711
Solubility (25°C)	3200ng/L
Sensitivity	7.18ng/L
Entrez	10664
Expiration Date	Please enquire.