

Product Datasheet

Human Tnfaip3-interacting Protein 1 ELISA Kit (orb1212838)

Description	Human Tnfaip3-interacting Protein 1, TNIP1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	TNIP1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	71,864
Uniprot ID	Q15025
Solubility (25°C)	3200ng/L
Sensitivity	7.52ng/L
Entrez	10318
Expiration Date	Please enquire.