

Product Datasheet

Human Thioredoxin-interacting Protein ELISA Kit (orb1212841)

Description	Human Thioredoxin-interacting Protein, TXNIP ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	TXNIP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	43,661
Uniprot ID	Q9H3M7
Solubility (25°C)	2400ng/L
Sensitivity	5.17ng/L
Entrez	10628
Expiration Date	Please enquire.