

Product Datasheet

Human Syndecan-2 ELISA Kit (orb1212852)

Description	Human Syndecan-2, SDC2 ELISA Kit
Reactivity	Human
Range	70-10000ng/L
Tested Applications	ELISA
Target	SDC2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,160
Uniprot ID	P34741
Solubility (25°C)	12800ng/L
Sensitivity	35.43ng/L
Entrez	6383
Expiration Date	Please enquire.