

Product Datasheet

Human Spondin-1 ELISA Kit (orb1212858)

Description	Human Spondin-1, SPON1 ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	SPON1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	90,973
Uniprot ID	Q9HCB6
Solubility (25°C)	24ng/mL
Sensitivity	0.028ng/mL
Entrez	10418
Expiration Date	Please enquire.