

Product Datasheet

Human Small Ubiquitin-related Modifier 2 ELISA Kit (orb1212863)

Description	Human Small Ubiquitin-related Modifier 2, SUMO2 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	SUMO2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	10,871
Uniprot ID	P61956
Solubility (25°C)	2400ng/L
Sensitivity	5.26ng/L
Entrez	6613
Expiration Date	Please enquire.