

Product Datasheet

Human Slam Family Member 7 ELISA Kit (orb1212865)

Description	Human Slam Family Member 7, SLAMF7 ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	SLAMF7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	37,421
Uniprot ID	Q9NQ25
Solubility (25°C)	9600ng/L
Sensitivity	14.56ng/L
Entrez	57823
Expiration Date	Please enquire.