

Product Datasheet

Human Septin-5 ELISA Kit (orb1212874)

Description	Human Septin-5, SEPT5 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	SEPTIN5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,777
Uniprot ID	Q99719
Solubility (25°C)	4800ng/L
Sensitivity	8.97ng/L
Entrez	5413
Expiration Date	Please enquire.