

Product Datasheet

Human Secretoglobin Family 1D Member 2 ELISA Kit (orb1212879)

Description	Human Secretoglobin Family 1D Member 2, SCGB1D2 ELISA Kit
Reactivity	Human
Range	20-4500ng/mL
Tested Applications	ELISA
Target	SCGB1D2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	9,925
Uniprot ID	095969
Solubility (25°C)	4800ng/mL
Sensitivity	9.32ng/mL
Entrez	10647
Expiration Date	Please enquire.