

Product Datasheet

Human Sal-like Protein 4 ELISA Kit (orb1212886)

Description	Human Sal-like Protein 4, SALL4 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	SALL4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	112,231
Uniprot ID	Q9UJQ4
Solubility (25°C)	4800ng/L
Sensitivity	8.91ng/L
Entrez	57167
Expiration Date	Please enquire.