

Product Datasheet

Human Roundabout Homolog 4 ELISA Kit (orb1212889)

Description	Human Roundabout Homolog 4, ROBO4 ELISA Kit
Reactivity	Human
Range	20-7000ng/L
Tested Applications	ELISA
Target	ROBO4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	107,457
Uniprot ID	Q8WZ75
Solubility (25°C)	8000ng/L
Sensitivity	9.26ng/L
Entrez	54538
Expiration Date	Please enquire.