

Product Datasheet

Human Ribonuclease 7 ELISA Kit (orb1212892)

Description	Human Ribonuclease 7, RNASE7 ELISA Kit
Reactivity	Human
Range	0.5-200ng/ml
Tested Applications	ELISA
Target	RNASE7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17,419
Uniprot ID	Q9H1E1
Solubility (25°C)	240ng/ml
Sensitivity	0.26ng/ml
Entrez	84659
Expiration Date	Please enquire.