

Product Datasheet

Human Retinoic Acid Receptor Beta ELISA Kit (orb1212899)

Description	Human Retinoic Acid Receptor Beta, RARB ELISA Kit
Reactivity	Human
Range	5-1000ng/L
Tested Applications	ELISA
Target	RARB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,489
Uniprot ID	P10826
Solubility (25°C)	1280ng/L
Sensitivity	2.16ng/L
Entrez	5915
Expiration Date	Please enquire.