

Product Datasheet

Human Ras-related Protein Rab-10 ELISA Kit (orb1212908)

Description	Human Ras-related Protein Rab-10, RAB10 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	RAB10
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,541
Uniprot ID	P61026
Solubility (25°C)	4800ng/L
Sensitivity	8.92ng/L
Entrez	10890
Expiration Date	Please enquire.