

Product Datasheet

Human Protein-glutamine Gamma-glutamyltransferase E ELISA Kit (orb1212915)

Description	Human Protein-glutamine Gamma-glutamyltransferase E, TGM3 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	TGM3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	76,632
Uniprot ID	Q08188
Solubility (25°C)	4800ng/L
Sensitivity	9.93ng/L
Entrez	7053
Expiration Date	Please enquire.