

Product Datasheet

Human Protein Wnt-16 ELISA Kit (orb1212921)

Description	Human Protein Wnt-16, WNT16 ELISA Kit
Reactivity	Human
Range	20-3500ng/L
Tested Applications	ELISA
Target	WNT16
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,690
Uniprot ID	Q9UBV4
Solubility (25°C)	4000ng/L
Sensitivity	9.61ng/L
Entrez	51384
Expiration Date	Please enquire.