

Product Datasheet

Human Protein S100-A7A ELISA Kit (orb1212923)

Description	Human Protein S100-A7A, S100A7A ELISA Kit
Reactivity	Human
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	S100A7A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,305
Uniprot ID	Q86SG5
Solubility (25°C)	48ng/mL
Sensitivity	0.058ng/mL
Entrez	338324
Expiration Date	Please enquire.