

Product Datasheet

Human Protein S100-A4 ELISA Kit (orb1212924)

Description	Human Protein S100-A4, S100A4 ELISA Kit
Reactivity	Human
Range	0.5-100ng/mL
Tested Applications	ELISA
Target	S100A4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,729
Uniprot ID	P26447
Solubility (25°C)	128ng/mL
Sensitivity	0.215ng/mL
Entrez	6275
Expiration Date	Please enquire.