

Product Datasheet

Human Protein S100-A2 ELISA Kit (orb1212925)

Description	Human Protein S100-A2, S100A2 ELISA Kit
Reactivity	Human
Range	0.5-200ng/ml
Tested Applications	ELISA
Target	S100A2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	11,117
Uniprot ID	P29034
Solubility (25°C)	240ng/ml
Sensitivity	0.29ng/ml
Entrez	6273
Expiration Date	Please enquire.