

Product Datasheet

Human Protein Kinase C-binding Protein Nell1 ELISA Kit (orb1212929)

Description	Human Protein Kinase C-binding Protein Nell1, NELL1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	NELL1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	89,635
Uniprot ID	Q92832
Solubility (25°C)	4800ng/L
Sensitivity	9.71ng/L
Entrez	4745
Expiration Date	Please enquire.