

Product Datasheet

Human Prostasin ELISA Kit (orb1212936)

Description	Human Prostasin, PRSS8 ELISA Kit
Reactivity	Human
Range	0.5-300ng/mL
Tested Applications	ELISA
Target	PRSS8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,431
Uniprot ID	Q16651
Solubility (25°C)	320ng/mL
Sensitivity	0.28ng/mL
Entrez	5652
Expiration Date	Please enquire.