

Product Datasheet

Human Procollagen C-endopeptidase Enhancer 1 ELISA Kit (orb1212950)

Description	Human Procollagen C-endopeptidase Enhancer 1, PCOLCE ELISA Kit
Reactivity	Human
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	PCOLCE
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	47,972
Uniprot ID	Q15113
Solubility (25°C)	48ng/mL
Sensitivity	0.056ng/mL
Entrez	5118
Expiration Date	Please enquire.