

Product Datasheet

Human Plasminogen Activator Inhibitor 2 ELISA Kit (orb1212958)

Description	Human Plasminogen Activator Inhibitor 2, SERPINB2 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	SERPINB2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	46,596
Uniprot ID	P05120
Solubility (25°C)	2400ng/L
Sensitivity	5.27ng/L
Entrez	5055
Expiration Date	Please enquire.