

Product Datasheet

Human Peptidyl-prolyl Cis-trans Isomerase C ELISA Kit (orb1212969)

Description	Human Peptidyl-prolyl Cis-trans Isomerase C, PPIC ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	PPIC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,763
Uniprot ID	P45877
Solubility (25°C)	3200ng/L
Sensitivity	7.36ng/L
Entrez	5480
Expiration Date	Please enquire.