

Product Datasheet

Human Pantetheinase ELISA Kit (orb1212975)

Description	Human Pantetheinase, VNN1 ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	VNN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	57,012
Uniprot ID	095497
Solubility (25°C)	9600ng/L
Sensitivity	14.55ng/L
Entrez	8876
Expiration Date	Please enquire.