

Product Datasheet

Human Nucleotide-binding Oligomerization Domain-containing Protein 2 ELISA Kit (orb1212979)

Description	Human Nucleotide-binding Oligomerization Domain-containing Protein 2, NOD2 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	NOD2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	115,283
Uniprot ID	Q9HC29
Solubility (25°C)	3200ng/L
Sensitivity	7.18ng/L
Entrez	64127
Expiration Date	Please enquire.