

Product Datasheet

Human Nucleoside Diphosphate Kinase A ELISA Kit (orb1212981)

Description	Human Nucleoside Diphosphate Kinase A, NME1 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	NME1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17,149
Uniprot ID	P15531
Solubility (25°C)	2400ng/L
Sensitivity	4.28ng/L
Entrez	4830
Expiration Date	Please enquire.