

Product Datasheet

Human Nucleophosmin ELISA Kit (orb1212982)

Description	Human Nucleophosmin, NPM1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	NPM1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	32,575
Uniprot ID	P06748
Solubility (25°C)	3200ng/L
Sensitivity	7.28ng/L
Entrez	4869
Expiration Date	Please enquire.