

Product Datasheet

Human Nucleobindin-1 ELISA Kit (orb1212984)

Description	Human Nucleobindin-1, NUCB1 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	NUCB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	53,879
Uniprot ID	Q02818
Solubility (25°C)	2400ng/L
Sensitivity	4.27ng/L
Entrez	4924
Expiration Date	Please enquire.