

Product Datasheet

Human Nuclease-sensitive Element-binding Protein 1 ELISA Kit (orb1212985)

Description	Human Nuclease-sensitive Element-binding Protein 1, YBX1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	RREB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	181,420
Uniprot ID	Q92766
Solubility (25°C)	4800ng/L
Sensitivity	9.29ng/L
Entrez	6239
Expiration Date	Please enquire.