

Product Datasheet

Human Nuclear Inhibitor of Protein Phosphatase 1 ELISA Kit (orb1212986)

Description	Human Nuclear Inhibitor of Protein Phosphatase 1, PPP1R8 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	PPP1R8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	38,479
Uniprot ID	Q12972
Solubility (25°C)	4800ng/L
Sensitivity	9.19ng/L
Entrez	5511
Expiration Date	Please enquire.