

Product Datasheet

Human Ninjurin-1 ELISA Kit (orb1212992)

Description	Human Ninjurin-1, NINJ1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	NINJ1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,345
Uniprot ID	Q92982
Solubility (25°C)	3200ng/L
Sensitivity	7.32ng/L
Entrez	4814
Expiration Date	Please enquire.