

Product Datasheet

Human Neuropilin-2 ELISA Kit (orb1212997)

Description	Human Neuropilin-2, NRP2 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	NRP2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	104,831
Uniprot ID	O60462
Solubility (25°C)	48ng/ml
Sensitivity	0.045ng/ml
Entrez	8828
Expiration Date	Please enquire.