

## Product Datasheet

### Human Neuritin ELISA Kit (orb1213001)

|                            |                                        |
|----------------------------|----------------------------------------|
| <b>Description</b>         | Human Neuritin, NRN1 ELISA Kit         |
| <b>Reactivity</b>          | Human                                  |
| <b>Range</b>               | 20-4500ng/L                            |
| <b>Tested Applications</b> | ELISA                                  |
| <b>Target</b>              | NRN1                                   |
| <b>Storage</b>             | 2-8°C                                  |
| <b>Note</b>                | For research use only                  |
| <b>Application notes</b>   | Sandwich                               |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates |
| <b>MW</b>                  | 15,333                                 |
| <b>Uniprot ID</b>          | <b>Q9NPD7</b>                          |
| <b>Solubility (25°C)</b>   | 4800ng/L                               |
| <b>Sensitivity</b>         | 9.61ng/L                               |
| <b>Entrez</b>              | <b>51299</b>                           |
| <b>Expiration Date</b>     | Please enquire.                        |