

Product Datasheet

Human Nephronectin ELISA Kit (orb1213005)

Description	Human Nephronectin, NPNT ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	NPNT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	61,907
Uniprot ID	Q6UXI9
Solubility (25°C)	24ng/mL
Sensitivity	0.029ng/mL
Entrez	255743
Expiration Date	Please enquire.