

Product Datasheet

Human Nadh-ubiquinone Oxidoreductase Chain 6 ELISA Kit (orb1213011)

Description	Human Nadh-ubiquinone Oxidoreductase Chain 6, MT-ND6 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	MT-ND6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,622
Uniprot ID	P03923
Solubility (25°C)	3200ng/L
Sensitivity	7.36ng/L
Entrez	4541
Expiration Date	Please enquire.