

## Product Datasheet

### Human Mucin-5Ac ELISA Kit (orb1213020)

|                            |                                           |
|----------------------------|-------------------------------------------|
| <b>Description</b>         | Human Mucin-5Ac, MUC5AC<br>ELISA Kit      |
| <b>Reactivity</b>          | Human                                     |
| <b>Range</b>               | 10-2000ng/mL                              |
| <b>Tested Applications</b> | ELISA                                     |
| <b>Target</b>              | MUC5AC                                    |
| <b>Storage</b>             | 2-8°C                                     |
| <b>Note</b>                | For research use only                     |
| <b>Application notes</b>   | Sandwich                                  |
| <b>Sample Types</b>        | Serum, plasma, cell culture<br>supernates |
| <b>MW</b>                  | 585,570                                   |
| <b>Uniprot ID</b>          | <b>P98088</b>                             |
| <b>Solubility (25°C)</b>   | 2400ng/mL                                 |
| <b>Sensitivity</b>         | 4.32ng/mL                                 |
| <b>Entrez</b>              | <b>4586</b>                               |
| <b>Expiration Date</b>     | Please enquire.                           |