

Product Datasheet

Human Melanoma-associated Antigen D1 ELISA Kit (orb1213035)

Description	Human Melanoma-associated Antigen D1, MAGED1 ELISA Kit
Reactivity	Human
Range	20-4500ng/L
Tested Applications	ELISA
Target	MAGED1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	86,161
Uniprot ID	Q9Y5V3
Solubility (25°C)	4800ng/L
Sensitivity	8.92ng/L
Entrez	9500
Expiration Date	Please enquire.