

Product Datasheet

Human Matrix Metalloproteinase-26 ELISA Kit (orb1213040)

Description	Human Matrix Metalloproteinase-26, MMP26 ELISA Kit
Reactivity	Human
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	MMP26
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,708
Uniprot ID	Q9NRE1
Solubility (25°C)	96ng/mL
Sensitivity	0.148ng/mL
Entrez	56547
Expiration Date	Please enquire.