

Product Datasheet

Human Leucine-rich Glioma-inactivated Protein 1 ELISA Kit (orb1213068)

Description	Human Leucine-rich Glioma-inactivated Protein 1, LGI1 ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	LGI1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	63,818
Uniprot ID	095970
Solubility (25°C)	24ng/mL
Sensitivity	0.021ng/mL
Entrez	9211
Expiration Date	Please enquire.