

Product Datasheet

Human Left-Right Determination Factor 2 ELISA Kit (orb1213069)

Description	Human Left-Right Determination Factor 2, LEFTY2 ELISA Kit
Reactivity	Human
Range	3-900ng/L
Tested Applications	ELISA
Target	LEFTY2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,920
Uniprot ID	000292
Solubility (25°C)	960ng/L
Sensitivity	1.63ng/L
Entrez	7044
Expiration Date	Please enquire.