

Product Datasheet

Human Laminin B1 ELISA Kit (orb1213071)

Description	Human Laminin B1, LN-B1 ELISA Kit
Reactivity	Human
Range	5-1000ng/mL
Tested Applications	ELISA
Target	LMNB1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	66,408
Uniprot ID	P20700
Solubility (25°C)	1280ng/mL
Sensitivity	2.74ng/mL
Entrez	4001
Expiration Date	Please enquire.