

Product Datasheet

Human Laminin Subunit Alpha-4 ELISA Kit (orb1213073)

Description	Human Laminin Subunit Alpha-4, LAMA4 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	LAMA4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	202,524
Uniprot ID	Q16363
Solubility (25°C)	3200ng/L
Sensitivity	7.96ng/L
Entrez	3910
Expiration Date	Please enquire.