

Product Datasheet

Human Krueppel-like Factor 4 ELISA Kit (orb1213075)

Description	Human Krueppel-like Factor 4, KLF4 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	KLF4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	54,671
Uniprot ID	043474
Solubility (25°C)	2560ng/L
Sensitivity	5.31ng/L
Entrez	9314
Expiration Date	Please enquire.