

Product Datasheet

Human Interleukin-31 Receptor Subunit Alpha ELISA Kit (orb1213082)

Description	Human Interleukin-31 Receptor Subunit Alpha, IL31RA ELISA Kit
Reactivity	Human
Range	3-900ng/L
Tested Applications	ELISA
Target	IL31RA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	82,954
Uniprot ID	Q8NI17
Solubility (25°C)	960ng/L
Sensitivity	1.34ng/L
Entrez	133396
Expiration Date	Please enquire.