

Product Datasheet

Human Interleukin-27 Subunit Beta ELISA Kit (orb1213083)

Description	Human Interleukin-27 Subunit Beta, EBI3 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	EBI3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,396
Uniprot ID	Q14213
Solubility (25°C)	2400ng/L
Sensitivity	4.64ng/L
Entrez	10148
Expiration Date	Please enquire.