

Product Datasheet

Human Interleukin-17 Receptor B ELISA Kit (orb1213091)

Description	Human Interleukin-17 Receptor B, IL17RB ELISA Kit
Reactivity	Human
Range	0.3-90ng/L
Tested Applications	ELISA
Target	IL17RB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	55,885
Uniprot ID	Q9NRM6
Solubility (25°C)	96ng/L
Sensitivity	0.17ng/L
Entrez	55540
Expiration Date	Please enquire.