

Product Datasheet

Human Interleukin-17 Receptor A ELISA Kit (orb1213092)

Description	Human Interleukin-17 Receptor A, IL17RA ELISA Kit
Reactivity	Human
Range	30-9000ng/L
Tested Applications	ELISA
Target	IL17RA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	96,122
Uniprot ID	Q96F46
Solubility (25°C)	9600ng/L
Sensitivity	16.15ng/L
Entrez	23765
Expiration Date	Please enquire.