

Product Datasheet

Human Interleukin-1 Receptor Accessory Protein ELISA Kit (orb1213098)

Description	Human Interleukin-1 Receptor Accessory Protein, IL1RAP ELISA Kit
Reactivity	Human
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	IL1RAP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	65,418
Uniprot ID	Q9NPH3
Solubility (25°C)	25.6ng/mL
Sensitivity	0.033ng/mL
Entrez	3556
Expiration Date	Please enquire.